
Read instructions carefully before starting test

Reveal[®] 3-D

for Gliadin

INTRODUCTION AND INTENDED USE

Reveal[®] 3-D for Gliadin is uniquely designed with 3 lines of detection and can be used virtually anywhere to screen environmental swabs and rinses for the presence of significant levels of gliadin (gluten). The test's 3-D technology ensures greater reliability with screening than ever before.

LIMIT OF DETECTION

The test demonstrates reactivity to both allergenic fractions of the gluten molecule, specifically gliadin and glutenin.

Utilizing the environmental swabs supplied, levels of 1 µg/100 cm² of gliadin extract can be detected.

When analyzing rinses, gliadin is detectable at a level of 5 ppm (10 ppm gluten). The presence of cleaners or sanitizers can affect the limit of detection (LOD) in rinses.

CROSS REACTIVITY

In addition to common bread wheat, the antibody used detects rye and barley.

Reveal 3-D for Gliadin is specific for the above cereals and will not detect oats, rice or maize.

A panel of potential cross reactants was assessed. No cross reactants were detected.

TEST PERFORMANCE

Reveal 3-D for Gliadin has undergone rigorous validation to evaluate the specificity, sensitivity, robustness and intra- and inter-batch variability of the test method on rinses and environmental swabs.

SAMPLE COMPATIBILITY

Reveal 3-D for Gliadin is designed to detect gliadin on environmental surfaces and in rinses. Gliadin is a protein found in wheat that belongs to a group of proteins called prolamins. Other prolamins include secalin, found in rye and hordein, found in barley. Reveal for Gliadin also detects secalin and hordein.

The validation of certain food ingredients may be applicable for the Reveal 3-D test format. However, it is recommended that a more broadly inclusive test method, such as Alert[®] or Veratox[®] for Gliadin R5, be used for food ingredients. Existing commodity validations should be revalidated when suppliers or the manufacturing process has been changed. Please contact a Neogen representative or distributor for additional details.

MATERIALS PROVIDED

The pack contains the following:

1. An instruction leaflet
2. 2 sealed foil pouches, each containing 10 green Reveal 3-D for Gliadin devices
3. 1 bottle of swab wetting solution
4. 20 bottles containing 20 mL of Extraction Buffer
5. 20 individually packaged, sterile swabs with break-off tips

SAMPLING TECHNIQUE

1. Clean-in-place rinse water

As only a small sample of material is required for Reveal 3-D for Gliadin, it is important to test a representative portion of the liquid.

2. Environmental swabbing

The swabs supplied are intended to be used for the collection of environmental samples from which the presence of gliadin (gluten) can be tested. This method can be used to validate the adequacy of cleaning and/or to identify problem areas (e.g., ineffective removal of gliadin (gluten) from processing equipment).

3. Food samples

Contact a Neogen representative or distributor for additional details.

PRECAUTIONS

1. This test should only be used for environmental testing in an industrial food manufacturing/preparation or for labeling enforcement.
2. Do not use any part of the test beyond the expiry date.
3. Do not open the foil bag until just before use and ensure it is tightly sealed after removal of a device.
4. Always store the kit between 2–8°C (35–46°F). Avoid freezing.

TEST STORAGE/SAMPLE STABILITY

Store the Reveal 3-D for Gliadin kit between 2–8°C (35–46°F), and use within the expiry date stated on the outer label. Do not freeze. Extracted samples should be used within **3 hours** of extraction.

LIMITATIONS

A negative test from a surface or rinse cannot exclude the possibility the food or swabbed surface has gliadin (gluten) present since it may be distributed unevenly on the surface or in the rinse and may be below the limit of detection of the test.

In-house validation should be performed to ensure the swabbing procedure can detect residual gliadin (gluten) at the desired levels.


Reveal 3-D for Gliadin is qualitative and should only be used as a preliminary screen for gliadin (gluten) content.

RINSE SAMPLING

Remove the appropriate number of the following and allow to equilibrate at room temperature before use (**20–30 minutes** out of refrigerator):

1. Extraction Buffer bottle
2. Reveal 3-D for Gliadin device (in foil pouch)


1. Add the 0.25 mL (250 μ L) sample to a new extraction buffer bottle.
2. Secure the cap and shake for **1 minute**.


SWAB SAMPLING

Remove the appropriate number of the following and allow to equilibrate at room temperature before use (**20–30 minutes** out of refrigerator):


1. Bottle of swab wetting solution
2. Extraction Buffer bottle
3. Reveal 3-D for Gliadin device (in foil pouch)
4. Sterile swab

1. Estimate a swabbing area of approximately 10 cm x 10 cm. Alternatively, use the swab to collect samples of contamination from problem areas (e.g., of processing equipment).


2. Gather the sample with the swab, using one of the following methods.


For dry surfaces: Remove a sterile swab from the packaging and wet with swab wetting solution. Swab a 10 x 10 cm area by using a crosshatch technique revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.

For wet surfaces: Remove a sterile swab from the pack-


aging and swab a 10 x 10 cm area by using a crosshatch technique revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing. Do not moisten swab prior to use.

3. Insert the swab to an Extraction Buffer bottle and carefully break off the moistened end at the prescored mark so that it remains in the bottle.
4. Secure the cap and shake for **1 minute**.


SAMPLE TESTING

1. Remove the lid and fill it with the liquid from the bottle. Any froth should remain in the bottle.
2. Dip the head of the Reveal device into the liquid in the lid. Ensure the cavity is saturated with the liquid.
3. Leave the cavity saturated until liquid is seen running in the test window.
4. Place device on a flat surface and allow test to develop for **5 minutes**.

NOTE: It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device.


READING RESULTS

Read the result **5 minutes** after dipping. Observations after 6 minutes may be inaccurate due to overdevelopment of the device. If distinct lines are not visible at positions O and C, refer to numbers 3 or 4 below.


1. Negative result

No line at position T (test): Level of gliadin/gluten below detection limit. (See **Limitations** section).



2. Positive result

Any intensity of line at position T (test): Level of gliadin/gluten above detection limit.

3. High positive results

No line is visible at position O (overload) and a line is faintly visible or absent at position T: Sample is overloaded with gliadin/gluten.

4. Invalid results

If no line appears at position C (control), then the test may be invalid.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's Web site at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

Please visit www.neogen.com/Corporate/termsconditions.html for Neogen's full terms and conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, crustacea, soy, walnut

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