



Captair Smart

Customized Instrument Enclosures



Filtered enclosures fully adapted to the size of your instruments, offering total protection to the user



www.erlab.com

Adapted to any kind of instruments generating chemicals !



Eelctrochemical Analyser



Mass spectrometer



HPLC



Thermal Analyser



Gas chromatograph



Dispersion analyser



Ion chromatograph



Easy to install and use!



Front access and rear access
through sliding doors.



Entry ports for
cables and hoses



High safety with the Captair Smart flexible filtration column




- A** Control panel
- B** Optional HEPA H14 filter
- C₁** Main molecular filter
- D** Fan box
- C₂** Safety molecular filter

- ☒ Main and safety molecular filters
- ☒ Optional HEPA filter for chemicals powders
- ☒ Optional filter monitoring system (for solvents)
- ☒ Filtration efficiency acc. to NFX 15-211 and JG/T385



Adaptable and durable construction!



- ☒ Dimensions L x D x H according to your requests!
- ☒ Structure made of black anodized 2,5 mm aluminium!
- ☒ Transparent panels made of acrylic 8 mm thickness
- ☒ Optional worktop made of black phenolic resin 12,5 mm thickness



GreenFumehoods
an alternative to
ducted fume hoods



Smart Fume Hoods
for your routine
handlings



Captair Smart
Weighing Stations



Captair Smart
Storage cabinets



Captair Store
Storage cabinets



Captair Flow
Clean Air Enclosures



Captair Bio
PCR Workstations



Captair Pyramid
Portable Glovebox



Halo
Laboratory Air
Filtration system



Halo Sense
Lab Air Quality Sensor

Europe

erlab D.F.S S.A.S (France)

Parc d'Affaires des Portes
BP 403
27104 Val de Reuil Cedex - France
Tel. : +33(0)2 32 09 55 80
Fax. : +33(0)2 32 09 55 90
E-Mail : Ventas@erlab.net

erlab S.L (Spain)

Pol. Ind. Sur
Pasaje Newton 3A
08754 El Papiol - Barcelona - Spain
Tel. : +34 93 673 24 74
Fax. : +34 93 673 24 76
E-Mail : Ventas@erlab.net

North America

erlab, Inc. (USA)

388 Newburyport Turnpike
Rowley, MA 01969 - USA
Tel. : +1(978) 948-2216
Fax. : +1(978) 948-3354
E-Mail : CaptairSales@erlab.com

Southeast Asia

erlab SEA Sdn Bhd (Malaysia)

25 Jalan Firma
2/1 Kawasan Perindustrian - Tebrau
81100 Johor Bahru - Malaysia
Tel. : +60 (0)7 3 555 724
Fax. : +60 (0)7 3 552 810
E-Mail : erlab@erlab.com.my

www.erlab.com

Asia, Middle East, Africa, South America

Kunshan erlab D.F.S Co.,Ltd

No.118 Feng Qin Rd.
KETD Kunshan
Jiangsu Province - China
Tel. : +86 (0)512 5781 4085
Fax. : +86 (0)512 5781 4082
E-Mail : Sales.china@erlab.com.cn

Your  distributor